

# Cellulite and Thermography

# INTRODUCTORY

Introductory professional thermographic system for detecting the thermal aspects of cellulite and choosing the most suitable treatments.



The CELL-METER mini PRO has been conceived as an INTRODUCTORY-PROFESSIONAL system, simple to use, to provide a quick approach to the modern method of thermographic screening.

# CELL-METER mini PRO

The use of contact thermography to detect the thermal aspects of cellulite is by now an integral part of the different treatments adopted against this unsightly condition. Ongoing technological and scientific development has given rise to numerous extremely versatile and reliable professional systems, ranging from applications in the medical sector to specific uses in beauty centres.

CELL-METER mini PRO is characterised by the following features:

- ease of use
- easy-to-interpret thermograms (thermal signs of cellulite)
- low cost of introduction to the method
- extreme versatility.

Given its specific characteristics, CELL-METER mini PRO cannot replace the specialised professional thermography systems of the CELLUVISION – THERMO-CELL-TEST – CELL-METER PROFESSIONAL lines, which represent the next step up for those who have become acquainted with this method thanks to the CELL-METER mini PRO.

## GENERAL CHARACTERISTICS AND POSSIBLE USES

The CELL-METER mini PRO as a whole may be defined as a multi-purpose system whose functions include:

- detection of the thermal signs of cellulite
- interpretation and classification of the stage of cellulite
- measurement of the leg circumference, for follow-up comparisons
- easy recording of the thermal signs detected

## COMPOSITION

The CELL-METER mini PRO is composed of:

- a thermographic plate including a tape measure in centimetres
- plate size: 12 x 70 cm
- temperature range covered: 4.8°C
- six-colour display of temperatures in the range between 28.3 and 33.1°C
- ELC thermographic surface area 100 sqcm
- protective container with 4-colour images
  - size when open: 45 x 28 cm
  - size when closed: 15 x 28 cm

including pockets for holding the thermographic plate and thermographic data recording charts.

- data recording charts
- poster to guide interpretation incorporated in the case
- user instructions in two languages (available versions: Italian and English or Spanish and Portuguese – other languages available on request)

## POTENTIAL USER TARGETS

There are numerous potential users of the CELL-METER mini PRO:

- **MANUFACTURERS** who already offer their treatment systems to professionals and wish to complete them with a diagnostic-educational tool.
- **MEDICAL OR BEAUTY ASSOCIATIONS** that intend to propose a simple introductory cellulite diagnosis system to their members
- **BEAUTY SCHOOLS**, which thanks to this practical tool can provide their students with an introductory system to acquaint them with the different aspects of cellulite.



## PROFESSIONAL AND SPECIALIST SYSTEMS



**THERMO-CELL-TEST DGT and MAC**  
Specialist systems with high-resolution thermographic plates and instant or digital imaging system for recording thermographic images.



**CELL-METER PROFESSIONAL KIT**  
System intended for professional use, comprising three thermographic plates with different thermal gradients and a practical tool for measuring leg circumference, as well as charts for recording the stages of cellulite detected.



**THERMO-CELL-TEST PROFESSIONAL KIT2**  
Easy-to-use instrument aimed at aestheticians and beauty centres, featuring two thermographic plates with different thermal gradients.



**CELLUVISION PROFESSIONAL KIT**  
Personalisable system featuring two thermographic plates, aimed at a broad target market given its relatively low cost.