



© 1993-2012 Inthesys 2000 Sas di V.Marzorati & C.

INTHESYS 2000 SAS di V.Marzorati e C.  
Viale Don Minzoni 38  
20091 BRESSO (Milano) Italy

Phone & fax: + 39.02.6650.3890  
Mobile phone : + 39.333.2409.797  
Fax server : + 39.02.9998.9476  
E-mail: [info@inthesys.it](mailto:info@inthesys.it)  
Codice fiscale e partita Iva : 10944320158  
C.C.I.A.A. Milano 1420476-Cod.mecc. :MI351680

## **MESOTHERAPY NEEDLES SPECIFICATIONS 2012**

<b>Product:</b> Mesotherapy needles
<b>Cannula:</b> Stainless steel AISI 304, raw, bright. drawn
<b>Connection for siryng:</b> In polipropylene resin hot printed
<b>Conicity:</b> Luer Lock, rules UNI 001/67
<b>Glue:</b> cannula is glued with Loctite FMD11/1811 kiln-dried
<b>Colours:</b> Various, it depends on Gauge international standard
<b>Sterilization:</b> with Etilene oxide gas
<b>Packaging:</b> single blister made of medical paper, porous ness controlled 100 needles/box
<b>Validity:</b> 5 years in sealed packaging
<b>Sizes 27 Gauge:</b> outside diameter 0.40 mm inside diameter 0.20 mm lenght 4mm (COD: <b>AAL4</b> ) - 6 mm (COD : <b>AAL6</b> ) (tolerance +/- 10%)
<b>Sizes 30 Gauge:</b> outside diameter 0.30 mm inside diameter 0.15 mm lenght 4mm (COD: <b>AAL34</b> ) - 6 mm (COD : <b>AAL36</b> ) (tolerance +/- 10%)
<b>Sizes 32 Gauge:</b> outside diameter 0.26 mm inside diameter 0.11 mm lenght 4 (COD: <b>AM 324</b> )- 6 (COD <b>AM 326</b> )- 12.7 mm (COD <b>AM32G</b> ) (tolerance +/- 10%)
<b>Distributor:</b> Inthesys 2000 SAS – 20091 BRESSO (Milano) – <a href="http://www.inthesys.it">www.inthesys.it</a>
<b>Manufactured:</b> In Italy – EC
<b>Conformity:</b> Directive 93/42 ECC (brought about Dlgs. 46 del 24/02/97) attached I
<b>Provisions:</b> EN 556, EN 550, EN 1441, ISO 10993-5, 10, 11, EN 1633-1
<b>CE certification:</b> N. 0476 (CERMET) EC CERTIFICATE N. MED-22009 valid until 2nd March 2007
<b>Class:</b> Classe II A